


LEGEND LP Tube is the latest bonding buccal tube created by TOMY Inc. LEGEND LP Tubes feature a large entry with chamfered edges to facilitate wire insertion along with gingival and occlusal side notches for a secure grip during positioning. LEGEND LP Tubes also include the same successful characteristics found in TOMY's current buccal tubes such as a curved anatomical contour base design, a "V" indent in crown, colour marks for easy quadrant identification and a micro-etched, 80-gauge mesh contour base for superior bonding. Its low profile and lowered mesially positioned hook enhances patient comfort.






Advantages

- Large mesial entry
- Comfortable
- Simple manipulation
- "V" indent in tooth crown for easy placement

LEGEND LP

Low Profile Bonding Buccal Tubes with Anatomical Micro-Etched Base



Experts in Orthodontics

GC Orthodontics Europe GmbH

www.gcorthodontics.eu


LEGEND LP

Quadrant ID Mark
Eliminates inversion

Chamfered Entry Edges
For easier wire insertion



Laser Welded Hooks



Gripping Notches
Assist with accurate positioning


"V" indent in crown for easy placement



Micro-Etched, 80-Gauge Mesh Base Ensures consistent bonding

Anatomical Base
With a "V" indent in crown

	Torque	Offset	R/L	P/N	
				.018"x.025"	.022"x.028"
U 6	-14°	10°	R	40-1160-0106	40-2160-0106
			L	40-1260-0106	40-2260-0106
U 6	-14°	0°	R	40-1160-0107	40-2160-0107
			L	40-1260-0107	40-2260-0107
U 7	-14°	10°	R	40-1170-0106	40-2170-0106
			L	40-1270-0106	40-2270-0106
U 7	-14°	0°	R	40-1170-0107	40-2170-0107
			L	40-1270-0107	40-2270-0107
L 6	-20°	0°	R	40-1460-0106	40-2460-0106
			L	40-1360-0106	40-2360-0106
L 7	-10°	0°	R	40-1470-0106	40-2470-0106
			L	40-1370-0106	40-2370-0106



99-40E0-13LP REV01-09.2019



Since 1921
Towards Century of Health

GC Orthodontics Europe GmbH
Harkortstraße 2 · 58339 Breckerfeld · Germany
Phone: +49 2338 801-888 · info.gco.germany@gc.dental